

DRI[®] Oxycodone Assay

IVD For In Vitro Diagnostic Use

REF	10015632 (3 x 18 mL Kit)
	100248 (70 mL Kit)
	100249 (500 mL Kit)

Intended Use

The DRI Oxycodone Enzyme Immunoassay is intended for the qualitative and semi-quantitative determination of oxycodone in human urine.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary results are used.

Summary and Explanation of the Test

Oxycodone is a semi-synthetic opioid prescribed for pain management in patients with moderate to severe pain. It is similar to codeine and morphine in its analgesic properties but it is more potent than morphine and has higher dependence potential. The drug oxycodone is supplied as OxyContin[®] (Oxycodone HCl) or in combination with aspirin (Percodan[®]) or acetaminophen (Percocet[®]).¹ Drug abusers crush the pills into powder and snort them for faster effect which may result in a potentially fatal outcome. According to Drug Abuse Warning Network (DAWN), there has been a dramatic increase in oxycodone related deaths.^{2,3} Oxymorphone, noroxycodone and noroxymorphone are the only known metabolites of oxycodone.² The metabolite, oxymorphone, is a potent narcotic analgesic, while the other two metabolites are relatively inactive. From 33-61% of a single dose of oxycodone is excreted in urine within 24 hours as unconjugated oxycodone (13-19%), conjugated oxycodone (7-29%), and conjugated oxymorphone (13-14%).⁴

The DRI[®] Oxycodone Assay is supplied as a liquid ready-to-use homogeneous enzyme immunoassay. The assay uses specific antibodies that can detect oxycodone and oxymorphone without any significant cross-reactivity to other opiate compounds. The assay is based on competition between a drug labeled with glucose-6-phosphate dehydrogenase (G6PDH), and free drug from the urine sample for a fixed amount of specific antibody binding sites. In the absence of free drug from the sample, the specific antibody binds the drug labeled with G6PDH and causes a decrease in enzyme activity. This phenomenon creates a direct relationship between the drug concentration in urine and enzyme activity. The enzyme activity is determined spectrophotometrically at 340 nm by measuring the conversion of nicotinamide adenine dinucleotide (NAD) to NADH.

Reagents

Antibody/Substrate Reagent:

Contains mouse monoclonal anti-oxycodone derivative antibody, glucose-6-phosphate (G6P), and nicotinamide adenine dinucleotide (NAD) in Tris buffer with sodium azide as a preservative.

Enzyme Conjugate Reagent:

Contains oxycodone derivative labeled with glucose-6-phosphate dehydrogenase (G6PDH) in Tris buffer with sodium azide as a preservative.

Additional Materials Required (sold separately):

REF	Kit Description
1664	DRI Negative Calibrator, 10 mL
1388	DRI Negative Calibrator, 25 mL
100250	DRI Oxycodone Calibrator 100, 10 mL
100251	DRI Oxycodone Calibrator 300, 10 mL
100252	DRI Oxycodone Calibrator 500, 10 mL
100253	DRI Oxycodone Calibrator 1000, 10 mL
DOAT-2	MAS [®] DOA Total – Level 2
DOAT-3	MAS [®] DOA Total – Level 3
DOAT-4	MAS [®] DOA Total – Level 4
DOAT-5	MAS [®] DOA Total – Level 5

⚠️ Precautions and Warnings

This test is for in-vitro diagnostic use only. The reagents are harmful if swallowed.

DANGER: DRI Oxycodone Assay contains ≤0.2% bovine serum albumin (BSA) and ≤0.5% Drug-specific antibody (Mouse).

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

Reagents used in the assay components contain ≤0.09% sodium azide. Avoid contact with skin and mucous membranes. Flush affected areas with copious amounts of water. Get immediate medical attention for eyes, or if ingested. Sodium azide may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up. Clean exposed metal surfaces with 10% sodium hydroxide.

Do not use the reagents beyond their expiration dates.

Reagent Preparation and Storage

The reagents are ready-to-use. No additional reagent preparation is required. The reagents should be stored refrigerated (2-8°C). All assay components, opened or unopened, are stable until the expiration date indicated on their respective labels. Do not use the reagents beyond their expiration dates.

Specimen Collection and Handling

Collect urine specimens in plastic or glass containers. Testing of fresh urine specimens is suggested.

The Mandatory Guidelines for Federal Workplace Drug Testing Programs; Final Guidelines recommend that specimens that do not receive an initial test within 7 days of arrival in the laboratory should be placed into secure refrigeration units.

Samples within a pH range of 3 to 11 are suitable for testing with this assay.

An effort should be made to keep pipetted samples free of gross debris. Centrifuge highly turbid specimens before analysis. Adulteration of the urine sample may cause erroneous results. If adulteration is suspected, obtain another sample and forward both specimens to the laboratory for testing.

Handle all urine specimens as if they were potentially infectious.

Assay Procedure

Analyzers capable of maintaining a constant temperature, pipetting samples, mixing reagents, measuring enzymatic rates at 340 nm and timing the reaction accurately can be used to perform this assay.

Refer to the specific application instructions of each analyzer for chemistry parameters before performing the assay.

Quality Control and Calibration

Qualitative analysis

For qualitative analysis of samples, use either the Oxycodone 100 Calibrator, or the Oxycodone 300 calibrator, as a cutoff level.

Semi-quantitative analysis

For semiquantitative analysis, use all calibrators.

Good laboratory practice suggests the use of control specimens to ensure proper assay performance. Use controls near the cutoff calibrator to validate the calibration. Control results must fall within established ranges, as determined by laboratory procedures and guidelines. If results fall outside of established ranges, assay results are invalid. All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements. Each laboratory should establish its own control frequency.

Results and Expected Values

Qualitative

Either the 100 or the 300 calibrator can be used as a cutoff reference for distinguishing "positive" from "negative" samples. A sample that exhibits a change in absorbance (ΔA) value equal to or greater than the value obtained with the cutoff calibrator is considered positive. A sample that exhibits a change in absorbance (ΔA) value lower than the value obtained with the cutoff calibrator is considered negative.

Semi-quantitative

A rough estimate of drug concentration in the samples can be obtained by running a standard curve with all calibrators and quantitating samples off the standard curve. Sample results above the high calibrator should be diluted with negative urine and retested.

Limitations

1. A positive result from this assay indicates only the presence of oxycodone or oxymorphone and does not necessarily correlate with the extent of physiological and psychological effects.
2. Care should be taken when reporting concentration results since there are many factors, e.g., fluid intake and other biologic factors, that may influence a urine test result.
3. It is possible that substances other than those investigated in the specificity study may interfere with the test and cause false results.

Typical Performance Characteristics

Typical performance data results obtained on the Hitachi 717 analyzer are shown below.⁵ The results obtained in your laboratory may differ from these data.

Precision

The DRI Oxycodone Controls (75, 125, 225, and 375 ng/mL) and cutoff calibrators (100 and 300 ng/mL) were tested in qualitative (mA) and semi-quantitative (ng/mL) mode using a modified NCCLS protocol. Results presented below were generated by testing all samples in replicates of 6, twice per day for 10 days using packaged reagents and controls on the Hitachi 717.

Qualitative (mA/min)

Calibrator/Control	100 ng/mL Cutoff					
	Within-run Precision			Total Precision		
	Mean	SD	% CV	Mean	SD	% CV
n=120						
75 ng/mL	348	2.1	0.6	348	2.9	0.8
100 ng/mL	371	1.9	0.5	371	3.2	0.9
125 ng/mL	389	2.0	0.5	389	3.1	0.8

Qualitative (mA/min)

Calibrator/Control	300 ng/mL Cutoff					
	Within-run Precision			Total Precision		
	Mean	SD	% CV	Mean	SD	% CV
n=120						
225 ng/mL	429	2.2	0.5	429	3.7	0.9
300 ng/mL	458	2.4	0.5	458	4.1	0.9
375 ng/mL	479	2.4	0.5	479	3.8	0.8

Semi-quantitative (ng/mL)

Calibrator/Control	100 ng/mL Cutoff					
	Within-run Precision			Total Precision		
	Mean	SD	% CV	Mean	SD	% CV
n=120						
75 ng/mL	73	2.4	3.3	73	2.9	4.0
100 ng/mL	98	2.9	2.9	98	3.6	3.7
125 ng/mL	123	2.4	2.4	123	4.9	4.0
225 ng/mL	227	5.0	2.2	227	8.2	3.6
300 ng/mL	303	9.0	3.0	303	11.5	3.8
375 ng/mL	375	10.2	2.7	375	14.7	3.9

Cutoff Characterization

Oxycodone samples around the cutoff were prepared by the addition of oxycodone stock solution to negative urine. The samples were targeted at the following control concentrations, 75 ng/mL and 125 ng/mL ($\pm 25\%$ of the 100 ng/mL cutoff) and 225 ng/mL and 375 ng/mL ($\pm 25\%$ of the 300 ng/mL cutoff). The samples were assayed in replicates of 21. Cutoff characterization was deemed acceptable if the observed oxycodone concentration for 95% of the 21 replicates was appropriately greater or lesser than the cutoff calibrator concentration. For all 21 replicates, the 75 ng/mL and 225 ng/mL samples assayed correctly, as less than their respective cutoff calibrators 100% of the time. The 125 ng/mL and 375 ng/mL samples assayed as greater than their respective cutoff calibrators 100% of the time.

Sensitivity

The sensitivity of the assay using the negative calibrator is 4.9 ng/mL.

Accuracy

One hundred and forty-four samples were analyzed by the DRI Oxycodone Assay in both the qualitative and semi-quantitative modes and the results were compared to the RapidOne™ OxyTest and to GC/MS. As the RapidOne™ OxyTest is a qualitative method for the detection of oxycodone at 100 ng/mL, only qualitative results at the 100 ng/mL cutoff for the DRI Oxycodone Assay were compared.

One hundred and twenty urine samples were assayed with the DRI Oxycodone Assay at 100 and 300 ng/mL cutoff on the Thermo Scientific Indiko and the Hitachi 717 analyzers.

Qualitative

The overall concordance between the DRI Oxycodone Assay and the RapidOne™ Oxy Test was 91.7%. The 12 samples detected as positive by the RapidOne™ Oxy Test and negative by the DRI Oxycodone Assay were confirmed by GC/MS to have oxycodone concentrations ≤ 100 ng/mL. At the 100 ng/mL cutoff, the overall concordance between DRI Oxycodone Assay and the GC/MS was 97.2%. At the 300 ng/mL cutoff, the overall concordance between DRI Oxycodone Assay and GC/MS was also 97.2%.

100 ng/mL Cutoff			100 ng/mL Cutoff			300 ng/mL Cutoff					
RapidOne™ Oxy			GC/MS			GC/MS					
	+	-		+	-		+	-			
DRI	+	65	0	DRI	+	61	4 [†]	DRI	+	39	1 [†]
	-	12	67		-	0	79		-	3	101

[†] Oxycodone concentrations ranged from 55-81 ng/mL.

[‡] Sample rate is 24 units above the cutoff in DRI Oxycodone Assay.

Semi-quantitative

The same 144 samples were assayed in tandem using GC/MS and the semi-quantitative mode of the DRI Oxycodone Assay. The results were analyzed using a 100 ng/mL cutoff.

GC/MS			
	+	-	
DRI	+	59	1*
	-	0	84

*Oxycodone concentration was 55 ng/mL by GC/MS with a DRI value of 103 ng/mL.

Specificity

The cross-reactivity of oxycodone metabolites, oxymorphone, noroxymorphone and noroxycodone, was evaluated by adding known amounts of each metabolite to oxycodone free urine. As indicated by the results in the table below, oxymorphone exhibits 103% cross reactivity with oxycodone; noroxymorphone and noroxycodone show no evidence of significant cross-reactivity.

Compound	Concentration tested (ng/mL)	Recovery (ng/mL)	% Cross-reactivity
Oxycodone	300	300	100
Oxymorphone	300	308	103
Noroxymorphone	500,000	303.5	<0.1
Noroxycodone	50,000	41.5	<0.1

The potential cross-reactivity posed by drugs commonly coadministered with oxycodone was evaluated by adding each substance to oxycodone free urine at the concentration indicated. A drug was considered to cross-react if the observed oxycodone concentration exceeded 100 ng/mL, the lowest cutoff for the DRI Oxycodone Assay. As shown in the tables below, all of the pharmacologic compounds evaluated, including a number of the opiate compounds, exhibited no cross-reactivity at the concentrations listed.

Structurally related opiate compounds that tested negative at 100 ng/mL cutoff.

Compound	Concentrations (µg/mL)
6-Acetyl Morphine	50
Codeine	500
Dihydrocodeine	100
Heroin	300
Hydrocodone	75
Hydromorphone	30
Levorphanol	200
Morphine	350
Morphine-3-glucuronide	900
Naloxone	200
Naltrexone	500
Norcodeine	1,000
Normorphine	1,000

Structurally unrelated compounds that tested negative at 100 ng/mL cutoff.

Compound	Concentrations (µg/mL)
Acetaminophen	1,000
Acetylsalicylic acid	1,000
Amitriptyline	500
Amoxicillin	500
Amphetamine	2,000
Benzoyllecgonine	2,000
Caffeine	1,000
Carbamazepine	1,000
Chlorpromazine	2,000
Clomipramine	1,000
Cimetidine	1,000
Desipramine	1,000
Dextromethorphan	200
Doxepine	200
Ephedrine	2,000
Fentanyl	200
Fluoxetine	1,000
Fluphenazine	500
Ibuprofen	1,000
Imipramine	1,000
Maprotiline	1,000
Meperidine	1,000
Methadone	1,000
Metroniazole	2,000
Nalbuphine	1,000
Nortriptyline	500
Oxazepam	500
Phencyclidine	1,000
Phenobarbital	1,000
Ranitidine	3,000
Secobarbital	1,000
Talwin	500
Thebaine	20
Thioridazine	1,000
Tramadol	500

Interference

The potential interference of pH and endogenous physiologic substances on recovery of oxycodone using the DRI Oxycodone Assay was assessed by spiking known amounts of potentially interfering substances into the low (225 ng/mL) and high (375 ng/mL) controls for the 300 ng/mL cutoff. No interference was observed by the addition of the compounds upto the concentrations listed below.

Compound	Concentrations (mg/dL)
Acetone	1,000
Ascorbic Acid	1,500
Creatinine	500
Ethanol	1,000
Galactose	10
Glucose	3,000
Hemoglobin	300
Human Serum Albumin	500
Oxalic Acid	100
Riboflavin	7.5
Sodium chloride	1,000
Urea	2,000
pH	3-11

References

1. Anderson D.T., Fritz K.L., and Muto J.J. OxyContin®: The concept of a "Ghost Pill" and the Postmortem Tissue Distribution of Oxycodone in 36 Cases. *J. Anal. Toxicol.* 2002, 26: 448-459.
2. *Clinical & Forensic Toxicology News*, Oxycodone: Recognition and Pharmacogenomics. By Jannetto P.J. and Gock S.B. March 2003.
3. Cone E.J., et al, Oxycodone Involvement in Drug Abuse Deaths: A DAWN-Based Classification Scheme applied to an Oxycodone Postmortem Database Containing over 1000 Cases. *J. Anal. Toxicol.* 2003, 27: 57-67.
4. Oxycodone. In: Baselt R.C. and Cravey R.H. *Disposition of toxic drugs and chemicals in man*, 4th ed. Chemical Toxicology Institute, Foster City, California: 1995: 572-574.
5. Data on file at Microgenics, a part of Thermo Fisher Scientific.



Microgenics Corporation
46500 Kato Road
Fremont, CA 94538 USA
US Customer and
Technical Support:
1-800-232-3342



EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau, Germany
Tel: +49 (0) 851 886 89 0
Fax: +49 (0) 851 886 89 10



For insert updates go to:
www.thermoscientific.com/diagnostics

Other countries:

Please contact your local Thermo Fisher Scientific representative.